

Arm sling

Shoulder brace for immobilisation in defined position

Indications

- preoperative
- postoperative
- posttraumatic
- subcapital fracture of the humerus when adequate fixation with strapping is not possible.

Features

- Partially elastic bandage material
- Individually adjustable strap



Mode of Operation

The **Arm Sling** is a brace to immobilise the shoulder in defined position. Shoulder immobilisation braces are made of strong textile support material, with Velcro fasteners. An appropriate strap restricts shoulder mobility in a defined position or immobilises the shoulder.

Size chart and article numbers

Size	Lenght of Arm sling	Art.-No.
S	38 cm	37-6000
M	44 cm	37-6001
L	50 cm	37-6002

Measurement Chart

Noch nicht richtig übersetzt

