

## Dynamic spring-loaded Hallux Valgus orthosis

**Patents** DE 10 2008 049 854

### Indications

- conservative treatment of moderate cases of hallux valgus as well as preventive treatment
- pre-surgical treatment in order to prepare the soft tissue for surgery
- post-surgical treatment in order to stabilize the surgery result

### Features

- therapy in the physiological position of the big toe
- **pain reduction**
- individual adjustment of the tension spring without tools
- individual adjustment of the straps
- easy handling and high wearing comfort
- available in four colours

**Optionally** available with a pair of socks.



### Mode of Operation

The halluxsan® is based on the CDS®-principle and has been designed to treat the Hallux Valgus. The brace applies a dynamic low load prolonged stretch to stimulate growth in the contracted tissue.



### Size chart and article numbers

colours	Art.-No. left	Art.-No. right
green	940A-L	940A-R
black	940AS-L	940AS-R
pink	940AP-L	940AP-R
grey	940AG-L	940AG-R

EU-shoe size	Art.-No. Socks left and right
35 - 38	70001
39 - 42	70002
43 - 46	70003

# Denotation of Measurements

