

Manusan

Finger Brace

Patents EP 3 421 016

Indications

- Joint contractures due to
 - burns
 - stroke
 - spinal injury
 - paralysis
 - craniocerebral injury
 - surgery
 - arthrosis and chronic polyarthritis
- cut injury
- tendon suture
- Morbus Dupuytren
- to prevent recurring contractures following arthrolysis

Features

- Extension and flexion therapy
- Individually adjustable spring force
- Individually adjustable shell and strap system
- Easy handling and comfortable in use
- Adjustable hand plate
- No tool necessary for adjusting spring force
- Individual toolless adaptation of chains to finger length



Mode of operation

The Manusan is used to treat an extension deficit of the finger joints. Active use trains the appendicular finger muscles. Passive use stimulates growth of contracted tissue due to continuous dynamic tension.

Size chart and article numbers

Size	L/L	M/M
Shells	Aluminium Shells	
Circumference of wrist [cm]	22 - 30	18 - 24
Width Palm plate [cm]	10	9
Art.-No. left	955LL-L	955MM-L
Art.-No. right	955LL-R	955MM-R

Accessories

	Art.-No. left	Art.-No. right
Thumb element	956	

Meaning of indicated dimension on page 124

Other sizes available on request

Other sizes available on request

Meaning of indicated dimensions

