

## Lumbar Bandage

### Indications

- Lower back pain
- Moderately severe lumbar sciatica with slight failures in intervertebral disc protrusion/prolapse
- Moderately severe radicular lumbar syndrome, pseudo-radicular lumbar syndrome
- Spondylolithesis, grade I with lumbago
- Moderately severe lumbar deformity (relapsing treatment, complicated course)
- Condition after intervertebral disc surgery with small up to medium-sized prolapse
- Spondylarthrosis
- Sacroiliac joint irritation
- Lumbar muscle insufficiency

### Features

- Additional lateral double straps
- Semi-elastic bandaging material
- Lateral reinforcing stays
- Firm abdominal piece with micro hook and loop fastener
- Bridging rods (spring band steel)



### Mode of operation

The MKS® Pontsana® functions according to the bridging principle. The wide bandage with the additional lateral double straps straightens and hence relieves the lumbar vertebral column.



### Size chart and article numbers

Size	1	2	3	4	5	6	51	61
<b>Body circumference [cm]</b>	60 - 70	70 - 80	80 - 90	90 - 105	105 - 120	120 - 135	105 - 120	120 - 135
<b>Abdominal height [cm]</b>	18	18	18	18	18	18	18	18
<b>Back height [cm]</b>	32	32	34	34	38	38	34 (reduced)	34 (reduced)
<b>Art.-No.</b>	582-801-S	582-802-S	582-803-S	582-804-S	582-805-S	582-806-S	582-8051-S	582-8061-S

#### Reduced height

<b>Abdominal height [cm]</b>	14	14	14	14	14	14
<b>Back height [cm]</b>	26	26	26	26	26	26
<b>Art.-No.</b>	582-801N-S	582-802N-S	582-803N-S	582-804N-S	582-805N-S	582-806N-S

# Dimension sheet torso

For orders please copy and send filled in form to:  
 E-Mail: info@albrechtgmbh.com  
 Fax: 08051-961 29 30

Order       Quotation      Date: \_\_\_\_\_

Medical supply store: \_\_\_\_\_

Customer ID: \_\_\_\_\_

Orthopaedic technician: \_\_\_\_\_

Patient ID: \_\_\_\_\_

## Product name

- MKS® Body belt
- MKS® Body belt plus
- MKS® Vario
- MKS® Lumbago
- MKS® Lumbago plus
- MKS® Lumbago Vario
- MKS® Pontsana®
- MKS® Pontsana® plus
- MKS® Pontsana® mobil
- MKS® Pontsana® stabil
- MKS® Osteo
- MKS® Osteo plus
- MKS® Osteo flex
- MKS® Osteo flex plus
- MKS® Lumbo
- MKS® Lumbo forte
- MKS® Lumbo plus
- MKS® Lumbo forte plus
- MKS® Lumbo mobil
- MKS® Lumbo fenestra
- MKS® Extensa
- MKS® Extensa plus
- MKS® Thorecta
- MKS® Thorecta plus
- MKS® Core Compression Bandage
- other: \_\_\_\_\_

**Abdominal height**  
 Maximum front height of the textile body belt

**Front height**  
 Measured from the upper edge of the pubic bone to the sternum (3-4 below the collarbone)

**Body circumference**  
 Measured at the middle of the torso (just above the navel)

**Back height reclination**

Measured from the sacrum (beginning of gluteal fold) to circa 2-3 cm below cervical vertebra C7

**Back height**

Measured from the sacrum to circa 2-3 cm below the tips of the shoulder blades



## Technical implementation (to be filled in by albrecht)

- Size modification
- Modular adaptation

Art.-No.: \_\_\_\_\_

Approval: \_\_\_\_\_